POLYAMIN-TEST

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Teil I

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Directions

Part I

1. Remove the glass cylinder from the kit, rinse it with the water sample and fill it up to the calibration mark.
2. Add reagent C until the base is well covered (approx. 25 drops).
3. Add 2 drops of reagent A and 2 drops of reagent B.
4. Close the glass cylinder with the screw cap and shake vigorously.
5. Wait until the solution in the cylinder separates into 2 layers.
Lower layer colourless: no polyamines present
Lower layer blue: polyamines present

Part II (when the lower layer is blue)

1. Aspirate titration solution from the storage bottle into the graduated pipette until the lower level of the black plunger ring is aligned with the top mark.
2. Open the glass cylinder, insert the full pipette and slowly add the solution in small amounts by pressing the plunger down.
3. After each addition close the glass cylinder, shake vigorously and wait until the lower layer has separated out.
4. Add further amounts until the base of the glass cylinder changes colour from blue to pale pink. Titration is then complete.
5. Read off the polyamine content in mg/l (or ppm) from the scale of the pipette at the lower level of the black plunger ring.
6. If the contents of one pipette are not sufficient to effect the colour change, fill the pipette with fresh titration solution as indicated in 1 above and continue the titration in the same sample according to 2 - 5. Calculate the result by adding the total amount of titration solution used for each test.

Caution: Reagent A contains sulphuric acid.
On contact with the eyes or skin, rinse with large amounts of water.

Reagent C contains chlorocarbons.
Do not allow it to enter the sewerage system!
Collect after use for disposal as harmful waste.

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